




## PROFILE® - II and PROFILE® - IIA with Paper Chain-of-Custody Forms



Training and Certification Program  
Presented by:

**MEDTOX®**

## Intended Use

The PROFILE<sup>®</sup> -II Drugs of Abuse Test is a one-step immunochromatographic test for the rapid, qualitative detection of one or more of the following: Amphetamines, Barbiturates, Benzodiazepines, Cocaine, Methamphetamine/ 3,4 Methylenedioxymethamphetamines, Opiates, Oxycodone, Phencyclidine, Propoxyphene and THC (Cannabinoids) in human urine. It is not for over the counter sale. (Specific configurations listed below.) Some product configurations include a Lateral Flow Adulteration panel (LFAS).

## PROFILE<sup>®</sup> -II / PROFILE<sup>®</sup> -IIA Panel Configurations

### PROFILE<sup>®</sup> -IIA 5 Panel

Detects the major metabolites of these drugs: Cannabinoids (11-nor-9-carboxy- $\Delta^9$ -THC), Cocaine (Benzoylecgonine), Opiates (Morphine-Codeine), Amphetamines (d-Amphetamine), and Phencyclidine and LFAS.

### PROFILE<sup>®</sup> -IIA 7 Panel

Detects the major metabolites of these drugs: Cannabinoids (11-nor-9-carboxy- $\Delta^9$ -THC), Cocaine (Benzoylecgonine), Opiates (Morphine-Codeine), Amphetamines (d-Amphetamine), Phencyclidine, Barbiturates (Phenobarbital), and Benzodiazepines (nordiazepam) and LFAS.

### PROFILE<sup>®</sup> -IIA 7 Panel

Detects the major metabolites of these drugs: Cannabinoids (11-nor-9-carboxy- $\Delta^9$ -THC), Cocaine (Benzoylecgonine), Opiates (Morphine-Codeine), Amphetamines (d-Amphetamine), Phencyclidine, Methamphetamine (mAMP), Oxycodone, and LFAS.

### PROFILE<sup>®</sup> -IIA 10 Panel

Detects the major metabolites of these drugs: Cannabinoids (11-nor-9-carboxy- $\Delta^9$ -THC), Cocaine (Benzoylecgonine), Opiates (Morphine-Codeine), Amphetamines (d-Amphetamine), Phencyclidine, Barbiturates (Phenobarbital), and Benzodiazepines (nordiazepam), Propoxyphene, and Methamphetamine (mAMP), Oxycodone and LFAS.

## The test detects drug classes at the following cutoff concentrations:

AMP	Amphetamine (d-Amphetamine)	1000 ng/ml
BAR	Barbiturates (Butalbital)	200 ng/ml
BZO	Benzodiazepines (Nordiazepam)	300 ng/ml
COC	Cocaine (Benzoylecgonine)	300 ng/ml
MAMP	Methamphetamine (d-methamphetamine)	1000 ng/ml
MDMA	3,4 Methylenedioxymethamphetamine (Ecstasy)	1500 ng/ml
OPI	Opiates (Codeine/Morphine)	2000 ng/ml
OXY	Oxycodone	100 ng/ml
PCP	Phencyclidine	25 ng/ml
PPX	Propoxyphene	300 ng/ml
THC	Cannabinoids (11-nor-9-carboxy- $\Delta^9$ -THC)	50 ng/ml

**LFAS** - Test devices contain a membrane strip laminated with adulterant test pads, testing for the presence of **Oxidants** and **Nitrites**, as well as determining approximate values of **pH** and **Specific Gravity**.

The PROFILE<sup>®</sup> -II/PROFILE<sup>®</sup> -IIA drugs of abuse test provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or High Performance Liquid Chromatography (HPLC) or Liquid Chromatography/Mass Spectrometry (LC/MS/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained.

**Kit Contents:**

The PROFILE<sup>®</sup> -II Drugs of Abuse Test System kit contains 25 test system bags and one instructional Quick Start guide. Product inserts are available at <http://www.medtox.com/ProductTraining.aspx>.

**Test System Contents:**

1. One test device and pipette in a foil package.
  - A. The test strips each contain a membrane coated with drug conjugate and a pad coated with antibody dye complexes in a protein matrix.
  - B. LFAS test devices contain a membrane strip laminated with adulterant test pads for testing the presence of Oxidants and Nitrites, as well as determining approximate values of Specific Gravity and pH in human urine.
2. Collection cup with temperature strip attached.
3. Split specimen urine transfer tubes.
4. One color comparator chart (for LFAS test devices).
5. Specimen biohazard bag.
6. UPS shipping supplies (lab pack, Saturday stickers and window pouch). To print your own UPS airbills, go to <http://medtox.upsrow.com>.
7. Customized onsite chain-of-custody forms.
8. Package insert.

**Precautions:**

1. Urine specimens and all materials coming in contact with them should be handled and disposed of as if infectious and capable of transmitting infection. Never pipette by mouth and avoid contact with broken skin.
2. Avoid cross-contamination of urine samples by using a new urine specimen container for each urine sample.
3. The device should remain in its original sealed foil pouch until ready to use. If the pouch is damaged, do not use the test.
4. Test kits should be stored at 36-77°F (2-25°C) and should be refrigerated if temperature exceeds 77°F (25°C).
5. If devices have been stored refrigerated, bring to ambient temperature (18-25°C / 64-77°F) prior to opening foil pouch.
6. Do not use tests after the expiration date printed on the package label.
7. The drug screen portion of the device is for in-vitro diagnostic use only. The LFAS strip is for forensic/toxicology use only.
8. PROFILE<sup>®</sup> -II provides only a preliminary, qualitative, analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Mass Spectrometry is the preferred confirmatory method (i.e. GC/MS or LC/MS/MS).
9. Adulteration may produce inaccurate test results.
10. Use a new pipette for each urine sample.
11. PROFILE<sup>®</sup> -II test device is for in-vitro diagnostic use only.
12. Do not overfill the sample well. Using more than 100 µl of urine in each sample well can cause over-saturation and may produce inaccurate results.

**Specimen Collection Procedures:**

- Step 1 - Verify donor's identity. Donor must provide a photo I.D. or be positively identified by employer's representative.
- Step 2 - Begin Step 1 of chain-of-custody form by having the donor write in their name and social security number/ employee id number. Then have the donor complete the MRO section on copy 4 of the chain-of-custody form.
- Step 3 - Donor's outer garments and personal belongings are to be left outside collection area. Wallet may be retained.
- Step 4 - Collector inspects the contents of the donor's pockets.
- Step 5 - Instruct donor to wash his hands prior to the collection.
- Step 6 - Open a MEDTOX test system bag and split specimen collection kit in the presence of the donor.
- Step 7 - Instruct donor to provide at least 60 milliliters of urine into the calibrated collection container with temperature strip.
- Step 8 - Check the temperature of the sample within four minutes. Acceptable range is 90° - 100° F (32° - 38° C). Record results in Step 2 on the chain-of-custody form.
- Step 9 - Inspect the sample for unusual odor, presence of foreign objects or material, or other signs of tampering (e.g. blue dye in specimen, smell of bleach).
- Step 10 - In full view of donor; both specimen transfer tubes are opened.
- Step 11 - Pour at least 30 milliliters of urine into each specimen transfer tube and secure lids.
- Step 12 - Place a security seal over the top and down the side of each urine transfer tube. Have the donor initial each security seal. Collector dates each seal.
- Step 13 - Complete Step 1 of chain-of-custody form by having donor sign the consent statement.
- Step 14 - Collector completes Step 3 by dating, printing and signing their name.
- Step 15 - Place the sealed specimens in the back pocket of the biohazard bag. Place all copies of the chain-of-custody form, except the green donor copy, into the front pocket of the biohazard bag. **DO NOT SEAL THE BAG.**
- Step 16 - Give the donor their copy of the chain-of-custody form. The collector procedure is now complete and the donor may leave the collection facility.
- Step 17 - The specimens are now secure for temporary storage.



**PROFILE<sup>®</sup> -II Test Procedures**

There are up to three drug test strips and lateral flow adulteration strip (on LFAS test devices) on the PROFILE<sup>®</sup>-II / PROFILE<sup>®</sup>-IIA devices. Each drug strip has a control (CTRL), which determines that the proper amount of urine has been used. The control line must form to have a valid result for that drug strip.


- Step 1 - Remove chain-of-custody form and one specimen from the biohazard bag.
- Step 2 - Remove foil sealed PROFILE<sup>®</sup>-II test device. Check the expiration date, record it and the lot number in Step 2 on the chain-of-custody form.
- Step 3 - Complete Step 3 of the chain-of-custody form by writing in the tester's name and obtaining their signature. (Only required if the tester and the collector of the sample are not the same individual.)
- Step 4 - Open one foil sealed device for the test to be performed.
- Step 5 - Hold the 100 µl sample pipette by the upper bulb.
- Step 6 - Open the sealed specimen bottle. Lower the pipette stem into the urine sample. Squeeze the upper bulb then release it slowly. (This motion will draw 100 µl of urine into the stem.) The urine sample should reach the top of the stem, and overflow slightly into the middle bulb. If not, repeat this process.
- Step 7 - Dispense the 100 µl sample into the sample well by squeezing the upper bulb. This will empty the stem delivering the entire contents into the sample well. Excess urine will remain in the overflow bulb. Apply 100 µl of sample ONE time into EACH sample well.
- Step 8 - Read pH, Specific Gravity, and Nitrites as soon as color changes on LFAS strip. Read oxidant at 60 seconds.
- Step 9 - Read the test at 5 minutes after sample application. Test results read after 5 minutes may not be consistent with those at 5 minutes.
- Step 10 - Complete the chain-of-custody form by removing the appropriate negative/further testing sticker and placing it onto the results sticker box in Step 4. (If applicable, turn to Copy 2 and annotate the results in the boxes provided.)
- Step 11 - If the test result is **NEGATIVE**:
- Fax the top copy of the chain-of-custody form to the MEDTOX Fax Result Line at 1-888-295-0466. (No cover page is necessary.)
  - Send the MRO copy (pink) to the MRO, if one is designated.
  - Discard the specimen and file the chain-of-custody form for your records.
- If the test needs **FURTHER TESTING (NON-NEGATIVE)**, or results are **Abnormal**:
- Place the top copy of the chain-of-custody form and the remaining sealed sample into the biohazard bag and transport to MEDTOX overnight via UPS. Apply Saturday deliver sticker to shipments being picked up on Fridays.
  - Send Copy 4 (pink) to the MRO, if one is designated, otherwise send it to the employer.
  - Discard the remaining unsealed specimen and file the clinic copy (yellow).

## Interpreting PROFILE<sup>®</sup> -II Test Results

### WHEN IN DOUBT - SEND IT OUT!


#### INVALID TEST

The absence of a reddish-purple Control (CTRL) line indicates the test is invalid. The urine sample should be retested on a new device. If the second test is also invalid, send the urine sample to a reference laboratory for additional testing.




#### NEGATIVE TEST

The appearance of both a reddish-purple Control (CTRL) line and a specific drug line indicates a negative test result. The color intensities of the Control line and a specific drug line may not be equal; any reddish-purple line visible at 5 minutes indicates a negative test result. Line intensity will vary from test to test.



#### NON-NEGATIVE TEST

The appearance of both a reddish-purple Control (CTRL) line and the absence of a line next to a specific drug name at 5 minutes indicates a preliminary positive test result for that drug. Occasionally a white line (line lighter than the background of the strip) may appear next to a specific drug name. This indicates a preliminary positive test result from the drug. In this example the COC and BZO tests are valid preliminary positives.



#### LFAS Test:

Urine that produces an abnormal result on the LFAS adulteration strip should be sent to a reference laboratory for more definitive testing to determine if the urine may be dilute, substituted, invalid and/or adulterated.

## Chain-of-Custody Form, Seals Intact

On Site Screening Custody Form - **FAX THIS COPY TO 888-295-0466**



402 W. County Rd D  
St. Paul, MN 55112  
(651) 636-7466  
(800) 832-3244

Employer:  
SAMPLE ACCOUNT  
SAMPLE NAME  
STREET ADDRESS  
CITY, STATE, ZIP CODE

MRO:  
DR JOHN SMITH  
STREET ADDRESS

CITY, MN ZIP CODE  
PH 651-636-7466

Account #



PJ399677

**STEP 1** To be completed by  
**COLLECTOR / DONOR**

Donors Name  
Last Name (or other I.D.)

Donor Social  
Security  
Number

Donor  
Phone

First  
Name

Referring Physician / Company

**DONOR CONSENT** I certify that I provided my specimen to the collector, that the specimen container was sealed with a tamper evident seal, and that the information provided on this form and on the label affixed to the specimen bottle is correct. I authorize MEDTOX to release the results of these tests to my employer or prospective employer and/or authorized healthcare professionals.

Signature

Date

MO DAY YEAR

PJ399677

To be Completed by **COLLECTOR** - Indicate Reason For Test

☐ Pre-employment ☐ Random ☐ Reasonable Suspicion ☐ Return to Duty  
☐ Follow-up ☐ Post Accident ☐ Other (specify):

Specimen temperature must be read within 4 minutes of collection.

Kit Exp. Date

Specimen temperature within range: ☐ Yes, 90° - 100°F / 32 - 38° C.  
☐ No, Record specimen temperature here

Lot #

COLLECTION SITE Phone No.

Date: Month Day Year

Collection Time

☐ AM  
☐ PM

I, the collector, by signing below certify that the specimen identified on this form is the specimen given to me by the donor identified above and that it has been collected, labeled, and sealed in accordance with applicable requirements. I also certify that I performed the on-site test unless signed by the tester.  
The tester certifies that the specimen was received from temporary storage and tested with standard procedures.

REMARKS CONCERNING  
COLLECTION / TEST

COLLECTOR

TESTER

Printed Name

Printed Name

Signature

Signature

TEST RESULTS:

Place Results  
Sticker Here



**MEDTOX LABORATORIES CONFIRMATION REQUEST:** (Send the marked laboratory specimen to MEDTOX Laboratories if result is non-negative.)

TEST(S) ORDERED

[ ] 1206

Complete Step 5 ONLY if sending specimen to the lab for testing.

DATE Month Day Year	RELEASED BY PRINTED NAME / SIGNATURE Tester's Printed Name Tester's Signature	RECEIVED BY PRINTED NAME / SIGNATURE <b>COURIER</b>	PURPOSE OF CHANGE For Transport to MEDTOX
------------------------	---	--	--

To be completed by **MEDTOX**

DATE Month Day Year	RELEASED BY PRINTED NAME / SIGNATURE <b>COURIER</b>	RECEIVED BY PRINTED NAME / SIGNATURE Received by Printed Name Received by Signature	PURPOSE OF CHANGE For Accessioning at MEDTOX
------------------------	--	---	---

100771 (8/08)

**COPY 1: IF ADDITIONAL TESTING IS NECESSARY, SEND WITH LABORATORY SPECIMEN TO MEDTOX**

SPECIMEN ID NO.

SPECIMEN ID NO.

PJ399677



PJ399677



echain<sup>®</sup> U.S. Patent No. US6,376,251 and US6,653,139

Date (Mo, Day, Yr)  
Donor's Initials



Further Testing

PJ399677



Negative



## Chain-of-Custody Form, Negative

On Site Screening Custody Form - **FAX THIS COPY TO 888-295-0466**



402 W. County Rd D  
St. Paul, MN 55112  
(651) 636-7466  
(800) 832-3244



**PJ399673**

Employer:

SAMPLE ACCOUNT  
SAMPLE NAME  
STREET ADDRESS  
CITY, STATE, ZIP CODE

MRO:

DR JOHN SMITH  
STREET ADDRESS

CITY, MN ZIP CODE  
PH 651-636-7466

Account #



**STEP 1 To be completed by COLLECTOR / DONOR**

Donors Name  
Last Name (or other I.D.)

DOE

Donor Social  
Security  
Number

9 8 7 6 5 4 3 2 1

Donor  
Phone

5 5 5 4 7 3 6 2 4 1

First  
Name

JANE

Referring Physician / Company

**DONOR CONSENT**

I certify that I provided my specimen to the collector, that the specimen container was sealed with a tamper evident seal, and that the information provided on this form and on the label affixed to the specimen bottle is correct. I authorize MEDTOX to release the results of these tests to my employer or prospective employer and/or authorized healthcare professionals.

Signature

Jane Doe

Date

2 / 24 / 09  
MO DAY YEAR

**PJ399673**

**STEP 2**

To be Completed by **COLLECTOR** - Indicate Reason For Test

Specimen temperature must be read within 4 minutes of collection.

☒ Pre-employment ☐ Random ☐ Reasonable Suspicion ☐ Return to Duty  
☐ Follow-up ☐ Post Accident ☐ Other (specify):

Kit Exp. Date

10 / 31 / 10

Specimen temperature ☒ Yes, 90° - 100°F / 32° - 38° C.

Lot #

12345

within range: ☐ No, Record specimen temperature here

Date: 02 / 24 / 2009

Collection Time

1030

☒ AM  
☐ PM

COLLECTION SITE Phone No.

8883242468

**STEP 3**

I, the collector, by signing below certify that the specimen identified on this form is the specimen given to me by the donor identified above and that it has been collected, labeled, and sealed in accordance with applicable requirements. I also certify that I performed the on-site test unless signed by the tester.  
The tester certifies that the specimen was received from temporary storage and tested with standard procedures.

**REMARKS CONCERNING COLLECTION / TEST**

COLLECTOR

TESTER

Printed Name JOHN SMITH

Printed Name

Signature John Smith

Signature

**STEP 4**

**TEST RESULTS:**



Negative 22222

**MEDTOX LABORATORIES CONFIRMATION REQUEST:** (Send the marked laboratory specimen to MEDTOX Laboratories if result is non-negative.)

TEST(S) ORDERED

I 1 1206

**STEP 5**

**Complete Step 5 ONLY if sending specimen to the lab for testing.**

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year	Tester's Printed Name Tester's Signature	<b>COURIER</b>	For Transport to MEDTOX

To be completed by **MEDTOX**

☐ SEAL INTACT

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year	<b>COURIER</b>	Received by Printed Name Received by Signature	For Accessioning at MEDTOX

100771 (8/08)

**COPY 1: IF ADDITIONAL TESTING IS NECESSARY, SEND WITH LABORATORY SPECIMEN TO MEDTOX**

eChain<sup>®</sup> U.S. Patent No: US 6,376,251 and US 6,653,139



## Chain-of-Custody Form, Further Testing

On Site Screening Custody Form - **FAX THIS COPY TO 888-295-0466**



402 W. County Rd D  
St. Paul, MN 55112  
(651) 636-7466  
(800) 832-3244



**PJ399671**

Employer:

SAMPLE ACCOUNT  
SAMPLE NAME  
STREET ADDRESS  
CITY, STATE, ZIP CODE

MRO:

DR JOHN SMITH  
STREET ADDRESS

CITY, MN ZIP CODE  
PH 651-636-7466

**STEP 1 To be completed by COLLECTOR / DONOR**

Donors Name  
Last Name (or other I.D.)

DOE

Donor Social  
Security  
Number

1 2 3 4 5 6 7 8 9

Donor  
Phone

5 5 5 6 9 8 6 1 3 9

First  
Name

JOHN

Referring Physician / Company

**DONOR CONSENT**

I certify that I provided my specimen to the collector, that the specimen container was sealed with a tamper evident seal, and that the information provided on this form and on the label affixed to the specimen bottle is correct. I authorize MEDTOX to release the results of these tests to my employer or prospective employer and/or authorized healthcare professionals.

Signature *John Doe*

Date 02/24/09  
MO DAY YEAR

Account #



**PJ399671**

**STEP 2**

To be Completed by **COLLECTOR** - Indicate Reason For Test

☒ Pre-employment ☐ Random ☐ Reasonable Suspicion ☐ Return to Duty  
☐ Follow-up ☐ Post Accident ☐ Other (specify):

Specimen temperature must be read within 4 minutes of collection.

Kit Exp. Date 10/31/10  
Lot # 12345

Specimen temperature ☒ Yes, 90° - 100°F / 32 - 38° C.  
within range: ☐ No, Record specimen temperature here

**STEP 3**

Date: 02/24/09

Collection Time 1030 AM

COLLECTION SITE Phone No.

8883242468

I, the collector, by signing below certify that the specimen identified on this form is the specimen given to me by the donor identified above and that it has been collected, labeled, and sealed in accordance with applicable requirements. I also certify that I performed the on-site test unless signed by the tester. The tester certifies that the specimen was received from temporary storage and tested with standard procedures.

REMARKS CONCERNING  
COLLECTION / TEST

COLLECTOR	TESTER
Printed Name JOHN SMITH	Printed Name
Signature <i>John Smith</i>	Signature

**STEP 4**

TEST RESULTS:



Further Testing 99999

**MEDTOX LABORATORIES CONFIRMATION REQUEST:** (Send the marked laboratory specimen to MEDTOX Laboratories if result is non-negative.)

TEST(S) ORDERED

[ ] 1206

**STEP 5**

Complete Step 5 ONLY if sending specimen to the lab for testing.

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year 02/24/09	Tester's Printed Name JOHN SMITH Tester's Signature <i>John Smith</i>	<b>COURIER</b>	For Transport to MEDTOX

To be completed by MEDTOX

☐ SEAL INTACT

**STEP 6**

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year	<b>COURIER</b>	Received by Printed Name Received by Signature	For Accessioning at MEDTOX

100771 (8/08)

**COPY 1: IF ADDITIONAL TESTING IS NECESSARY, SEND WITH LABORATORY SPECIMEN TO MEDTOX**

eChain<sup>®</sup> U.S. Patent No. US-6,376,251 and US-6,653,139

## Chain-of-Custody Form, MRO

On Site Screening Custody Form



402 W. County Rd D  
St. Paul, MN 55112  
(651) 636-7466  
(800) 832-3244



PJ399671

Employer:

SAMPLE ACCOUNT

SAMPLE NAME

STREET ADDRESS

CITY, STATE, ZIP CODE

MRO:

DR JOHN SMITH

STREET ADDRESS

CITY, MN ZIP CODE

PH 651-636-7466

Account #



**STEP 1** To be completed by  
**COLLECTOR / DONOR**

Donor Social  
Security  
Number

1 2 3 4 5 6 7 8 9

Donors Name  
Last Name (or other I.D.)

Donor  
Phone

5 5 5 6 9 8 6 1 3 9

DOE

First  
Name

JOHN

Referring Physician / Company

**DONOR CONSENT**

I certify that I provided my specimen to the collector, that the specimen container was sealed with a tamper evident seal, and that the information provided on this form and on the label affixed to the specimen bottle is correct. I authorize MEDTOX to release the results of these tests to my employer or prospective employer and/or authorized healthcare professionals.

Signature John Doe

Date 02/24/09  
MO DAY YEAR

**STEP 2**

To be Completed by **COLLECTOR** - Indicate Reason For Test

Specimen temperature must be read within 4 minutes of collection.

☒ Pre-employment ☐ Random ☐ Reasonable Suspicion ☐ Return to Duty  
☐ Follow-up ☐ Post Accident ☐ Other (specify):

Kit Exp. Date 10/31/10  
Lot # 12345

Specimen temperature ☒ Yes, 90° - 100°F / 32 - 38° C.  
within range: ☐ No, Record specimen temperature here

**STEP 3**

Date: 02/24/09

Collection Time 1030 ☒ AM ☐ PM

COLLECTION SITE Phone No.

8883242468

I, the collector, by signing below certify that the specimen identified on this form is the specimen given to me by the donor identified above and that it has been collected, labeled, and sealed in accordance with applicable requirements. I also certify that I performed the on-site test unless signed by the tester.  
The tester certifies that the specimen was received from temporary storage and tested with standard procedures.

REMARKS CONCERNING  
COLLECTION / TEST

COLLECTOR	TESTER
Printed Name <u>JOHN SMITH</u>	Printed Name
Signature <u>John Smith</u>	Signature

To be completed by Donor

Donor Name JOHN DOE

Donor Address 1234 FIFTH AVENUE

City ANYTOWN

State UT

Zip 54321

Home Phone # (555) 698-6139

Work Phone # (555) 456-1237

To be completed by the **MEDICAL REVIEW OFFICER**

I have reviewed the laboratory results for the specimen identified by this form in accordance with applicable requirements. My determination / verification is:

☐ Negative ☐ Positive ☐ Test Not Performed ☐ Test Cancelled

REMARKS:

(PRINT) Medical Review Officer's Name (First, MI, Last)

Signature of Medical Review Officer

Date (Mo., Day, Yr.)



Complete Step 5 ONLY if sending specimen to the lab for testing.

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year <u>2/24/09</u>	Tester's Printed Name <u>John Smith</u> Tester's Signature	<b>COURIER</b>	For Transport to MEDTOX

To be completed by **MEDTOX**

☐ SEAL INTACT

DATE	RELEASED BY PRINTED NAME / SIGNATURE	RECEIVED BY PRINTED NAME / SIGNATURE	PURPOSE OF CHANGE
Month Day Year <u>2/24/09</u>	<b>COURIER</b>	Received by Printed Name Received by Signature	For Accessioning at MEDTOX

100771 (8/08)

COPY 4: SEND TO MRO (IF LISTED ABOVE) OTHERWISE EMPLOYER

eChain<sup>®</sup> U.S. Patent No. US 6,376,251 and US 6,653,139

**Limitations of Procedure**





- The PROFILE<sup>®</sup>-II Drugs of Abuse Test System is only for use with unadulterated human urine samples. Urine samples which are either extremely acidic (below pH 4.0) or basic (above pH 9.0) may produce erroneous results.
- A non-negative (presumptive positive) test result for a specific drug indicates that the sample may contain drug/drug metabolite near or above the cutoff level. It does not indicate the level of intoxication or the specific concentration of drug in the urine sample. Non-negative (presumptive positive) samples should be sent to a reference laboratory for more definitive testing.
- Test results interpreted after 5 minutes may not be reliable.
- The PROFILE<sup>®</sup>-II Drugs of Abuse Test System was not evaluated in point-of-collection settings.
- There is a possibility that other substances and/or factors, e.g. technical or procedural errors, may interfere with the test and cause false results.

**You have now completed the PROFILE<sup>®</sup>-II training program. To achieve certification as a tester with this device, you must successfully complete the following ten question certification quiz with a score of 80% or higher.**

**If you have any questions and would like to speak to a service representative, please call 888-324-2468 or email us at [echain@medtox.com](mailto:echain@medtox.com).**



# PROFILE® - II and PROFILE® - IIA with Paper Chain-of-Custody Forms




Certification Quiz  
Presented by:

**MEDTOX®**

## PROFILE<sup>®</sup> -II CERTIFICATION QUIZ


After reading each question completely, choose the best answer. Record your answers on the answer sheet found at the end of the quiz. To receive a certification of achievement, fax your answer sheet to MEDTOX at 866-398-3783, or email to [echain@medtox.com](mailto:echain@medtox.com).

1. Depending on configuration, PROFILE<sup>®</sup> -II tests a single urine sample for:
  - A. amphetamines, cocaine, opiates, PCP, THC
  - B. amphetamines, barbiturates, benzodiazepines, cocaine, opiates, PCP, THC
  - C. amphetamines, barbiturates, benzodiazepines, cocaine, methamphetamine, opiates, Oxycodone, PCP, propoxyphene, THC
  - D. All of the above
2. The entire test is invalid if:
  - A. No control (CTRL) line forms in the window
  - B. The urine sample smells bad
  - C. Control (CTRL) line only forms
  - D. A and C only
3. PROFILE<sup>®</sup> -II devices should be stored at:
  - A. 25° -57° F
  - B. 36° -77° F
  - C. 25° -77° F
  - D. 35° -88° F
  - E. 15° -57° F
4. PROFILE<sup>®</sup> -II provides a:
  - A. Confirmed analytical result
  - B. Quantitative analytical result
  - C. Preliminary test result
  - D. Qualitative test result
  - E. C and D only
5. Interpret the following test:
  - A. Negative for all five drugs
  - B. Valid test
  - C. Invalid test
  - D. Non-negative for THC
  - E. A and B only




6. Interpret the following test:

- A. Non-Negative for amphetamines
- B. Negative for all seven drugs
- C. Valid test
- D. Invalid test
- E. A and C only




7. Interpret the following test:

- A. Non-Negative for THC, PCP
- B. Negative for THC, PCP
- C. Negative for all nine drugs
- D. Non-negative for OPI, AMP, COC
- E. Invalid test



8. Interpret the following test:

- A. Non-Negative for THC, COC
- B. Non-Negative for OPI, AMP, PCP
- C. Negative for OPI only
- D. Negative for THC, COC
- E. Invalid test





9. The control (CTRL) line:
- A. Should exhibit a reddish-purple line
  - B. Should have no line appear
  - C. Is used as a comparative indicator
  - D. B and C only
10. PROFILE<sup>®</sup> -II test results may be safely interpreted:
- A. Immediately after the test is completed
  - B. As soon as the control (CTRL) line appears in the window
  - C. At 5 minutes
  - D. Within one hour after the test is completed
11. The test kit may be used after the expiration date:
- A. If the kit has not been opened
  - B. Never
  - C. Within one month after the expiration date
  - D. Only if the kit has been frozen
12. Each test result should be:
- A. Separately and independently compared to the control line
  - B. Separately and independently compared to the negative control line
  - C. Separately and independently evaluated for the presence or absence of a line
  - D. Compared to adjacent test lines
13. The following information is required on the chain-of-custody form:
- A. Donor identification number
  - B. Date of collection
  - C. Specimen temperature in or out of range
  - D. Name of collector and tester
  - E. All of the above
14. When the donor presents the urine sample to the collector:
- A. Donor must leave the facility prior to dividing the sample into the urine transfer tubes
  - B. The collector will already have opened specimen transfer tubes and the collector will be ready to split the sample
  - C. Collector should check the temperature of the sample and annotate it on the chain-of-custody form in Step 2
  - D. Collector should open both urine transfer tubes in full view of the donor and pour at least 30 ml into each urine transfer tube from the original collection container
  - E. Both C and D
15. The following statement is not true:
- A. Donor is instructed to wash hands prior to collection
  - B. Donor is instructed to provide at least 60 ml of urine
  - C. Donor initials and dates the security seals that have been affixed to the urine transfer tubes
  - D. Donor is instructed that he cannot wear his coat inside the collection area
  - E. Photo I.D. is not required if positive identification is provided by the employer's representative

**PROFILE<sup>®</sup> -II****Certification Quiz Answer Sheet**

Enter your answers in the spaces below, fill in your name and other requested information, and fax your answer sheet to MEDTOX at 866-398-3783, or email to [echain@medtox.com](mailto:echain@medtox.com). A Certification of Achievement will be faxed to you if you receive a score of 80% or higher.

Question 1: \_\_\_\_\_

Question 6: \_\_\_\_\_

Question 11: \_\_\_\_\_

Question 2: \_\_\_\_\_

Question 7: \_\_\_\_\_

Question 12: \_\_\_\_\_

Question 3: \_\_\_\_\_

Question 8: \_\_\_\_\_

Question 13: \_\_\_\_\_

Question 4: \_\_\_\_\_

Question 9: \_\_\_\_\_

Question 14: \_\_\_\_\_

Question 5: \_\_\_\_\_

Question 10: \_\_\_\_\_

Question 15: \_\_\_\_\_

Name (please print): \_\_\_\_\_

Company Name: \_\_\_\_\_

Address: \_\_\_\_\_

City, State, Zip : \_\_\_\_\_

Telephone: \_\_\_\_\_

Fax: \_\_\_\_\_

Email: \_\_\_\_\_

Or Mail to: MEDTOX Laboratories  
Attn: PROFILE<sup>®</sup> Certification  
402 West County Road D  
Saint Paul, MN 55112

Or email answers to: [echain@medtox.com](mailto:echain@medtox.com)

Please call 888-324-2468 with any questions, or email [echain@medtox.com](mailto:echain@medtox.com).

8.2009

© Copyright 2009, MEDTOX Scientific, Inc.

MEDTOX<sup>®</sup>, PROFILE<sup>®</sup> -II and PROFILE<sup>®</sup> -IIA are registered trademarks.