



EZ-SCREEN[®] (LO) Cup

DRUGS OF ABUSE TESTING

DRUG TESTING MADE EZ

MEDTOX[®] makes testing for drugs of abuse extremely easy — simply collect a sample and read the screening results. Our EZ-SCREEN Cup provides fast and accurate screening results in minutes.

EZ-SCREEN (LO) Cup Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485:2003
- Accurate screening results in minutes
- On-board adulteration strip available on some products
- Panel options that can test for up to 11 drugs — including Oxycodone and Methamphetamine/Ecstasy
- Technical support services available for all MEDTOX Diagnostics, Inc. products

EZ-SCREEN (LO) Screening Thresholds

Amphetamine (AM)	-	300 ng/ml
Barbiturates (BA)	-	200 ng/ml
Benzodiazepine (BZ)	-	200 ng/ml
Buprenorphine (BU)	-	10 ng/ml
Cocaine (CO)	-	100 ng/ml
Methamphetamine (MA)	-	1000 ng/ml
Methadone (MT)	-	200 ng/ml
Opiates (OP)	-	100 ng/ml
Oxycodone (OX)	-	100 ng/ml
Phencyclidine (PC)	-	25 ng/ml
Cannabinoids (TH)	-	40 ng/ml

EZ-SCREEN Cup Test Panel Options

Part No.	Product Description	Panel Configurations
605104	EZ-SCREEN (LO) Cup 5 Panel	AM/BZ/CO/OP/TH
605101	EZ-SCREEN (LO) Cup 6 Panel	AM/BZ/CO/MA/OP/TH
605103	EZ-SCREEN (LO) Cup 7 Panel	AM/BZ/CO/MA/OP/OX/TH + AS
605100	EZ-SCREEN (LO) Cup 9 Panel	AM/BZ/CO/MA/MT/OP/OX/PC/TH
605102	EZ-SCREEN (LO) Cup 11 Panel	AM/BA/BU/BZ/CO/MA/MT/OP/OX/PC/TH + AS

MEDTOX controls must be used with MEDTOX devices.