



## EZ-SCREEN<sup>®</sup> (HI) & PROFILE<sup>®</sup>-IIA DRUGS OF ABUSE TESTING

### FLEXIBLE AND QUICK RESULTS.

With MEDTOX<sup>®</sup> drug testing devices, you have a choice of a cup or cassette. EZ-SCREEN (HI) cup or our PROFILE-IIA cassette provide fast and accurate screening results in minutes.

#### EZ-SCREEN or PROFILE-IIA Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485:2003
- Accurate screening results in minutes
- On-board adulteration strip on all EZ-SCREEN (HI) and PROFILE-IIA devices
- Panel options that can test for up to 11 drugs — including Oxycodone and Methamphetamine/Ecstasy
- Technical support services available for all MEDTOX Diagnostics, Inc. products

#### EZ-SCREEN and PROFILE-IIA Screening Thresholds

Amphetamine (AM/AMP)	- 1000 ng/ml
Barbiturates (BA/BAR)	- 200 ng/ml
Benzodiazepine (BZ/BZO)	- 300 ng/ml
Cocaine (CO/COC)	- 300 ng/ml
Methamphetamine (MA/MAMP)	- 1000 ng/ml
Methadone (MT/MTD)	- 300 ng/ml
Opiates (OP/OPI)	- 2000 ng/ml
Oxycodone (OX/OXY)	- 100 ng/ml
Phencyclidine (PC/PCP)	- 25 ng/ml
Propoxyphene (PP/PPX)	- 300 ng/ml
Cannabinoids (TH/THC)	- 50 ng/ml

#### PROFILE-IIA Test Panel Options

Part No.	Product Description	Panel Configurations
601599	PROFILE-IIA 5 Panel	AMP/COC/OPI2/PCP/THC + LFAS
601607	PROFILE-II 7 Panel (BAR/BZO)	AMP/BAR/BZO/COC/OPI2/PCP/THC
601673	PROFILE-IIA 7 Panel (MAMP/OXY)	AMP/COC/MAMP/OPI2/OXY/PCP/THC + LFAS
601677	PROFILE-IIA 10 Panel	AMP/BAR/BZO/COC/MAMP/OPI2/OXY/PCP/PPX/THC + LFAS

#### EZ-SCREEN Cup Test Panel Options

Part No.	Product Description	Panel Configurations
605200	EZ-SCREEN (HI) 5 Panel	AM/CO/OP/PC/TH + AS
605201	EZ-SCREEN (HI) 6 Panel	AM/CO/MA/OP/PC/TH + AS
605202	EZ-SCREEN (HI) 7 Panel	AM/CO/MA/OP/OX/PC/TH + AS
605203	EZ-SCREEN (HI) 8 Panel	AM/BA/BZ/CO/MA/OP/PC/TH + AS
605204	EZ-SCREEN (HI) 10 Panel	AM/BA/BZ/CO/MA/MT/OP/OX/PC/TH + AS
605205	EZ-SCREEN (HI) 11 Panel	AM/BA/BZ/CO/MA/MT/OP/OX/PC/PP/TH + AS

*MEDTOX controls must be used with MEDTOX devices.*