



PROFILE-II ER Screening Thresholds

Amphetamine (AMP)	-	1000 ng/ml
Barbiturates (BAR)	-	200 ng/ml
Benzodiazepine (BZO)	-	300 ng/ml
Cocaine (COC)	-	300 ng/ml
Methamphetamine (MAMP)	-	1000 ng/ml
3,4-Methylenedioxy-N-		
Methylamphetamine (MDMA)	-	1500 ng/ml
Methadone (MTD)	-	300 ng/ml
Opiates (OPI3)	-	300 ng/ml
Oxycodone (OXY)	-	100 ng/ml
Phencyclidine (PCP)	-	25 ng/ml
Propoxyphene (PPX)	-	300 ng/ml
Tricyclic Antidepressant (TCA)	-	300 ng/ml
Cannabinoids (THC)	-	50 ng/ml

PROFILE-II ER[®]

DRUGS OF ABUSE TESTING

PERFORMANCE. EFFICIENCY. VALUE.

MEDTOX[®] PROFILE-II ER drug testing devices are accurate, easy to use and cost effective. Simply add the urine specimen to one of our proven cassette devices and get accurate screening results in minutes.

PROFILE-II ER Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485:2003
- Accurate screening results in minutes
- Panel options that can test for up to 12 drugs — including Oxycodone and Methamphetamine/Ecstasy
- Technical support services available for all MEDTOX Diagnostics, Inc. products

Ordering Information for PROFILE-II ER Devices

Part No.	Product Description	Panel Configurations
601591	PROFILE-II ER 9 Panel	AMP, BAR, BZO, COC, MTD, OPI, PCP, TCA, THC
601636	PROFILE-II ER 10 Panel	AMP, BAR, BZO, COC, MAMP (MDMA), MTD, OPI3, PPX, TCA, THC
601672	PROFILE-II ER 11 Panel	AMP, BAR, BZO, COC, MAMP (MDMA), MTD, OPI3, PCP, PPX, TCA, THC
601674	PROFILE-II ER 12 Panel	AMP, BAR, BZO, COC, MAMP (MDMA), MTD, OPI, OXY, PCP, PPX, TCA, THC <i>This device includes Oxycodone</i>

MEDTOX controls must be used with MEDTOX devices.