



VERDICT[®]-II

DRUGS OF ABUSE TESTING

QUICK, ACCURATE RESULTS

Rely on VERDICT-II cassette devices for reliable drugs of abuse testing. These devices deliver accurate screening results within minutes.

VERDICT-II Screening Thresholds

Amphetamine (AMP)	- 1000 ng/ml
Barbiturates (BAR)	- 200 ng/ml
Benzodiazepine (BZO)	- 300 ng/ml
Buprenorphine (BUP)	- 10 ng/ml
Cannabinoids (THC)	- 50 ng/ml
Cocaine (COC)	- 300 ng/ml
Opiates (OPI3)	- 300 ng/ml
Oxycodone (OXY)	- 100 ng/ml
Phencyclidine (PCP)	- 25 ng/ml
Tricyclic Antidepressant (TCA)	- 300 ng/ml

VERDICT-II Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485
- Accurate screening results in minutes
- Panel options include up to 7 tests in a single device
- Technical support services available for all MEDTOX Diagnostics, Inc. products

VERDICT-II — Traditional Drug Cut-Off Level Cassette with Pipette

Part No.	Product Description	Panel Configurations
601629	VERDICT-II	TCA
601675	VERDICT-II	Oxycodone
601682	VERDICT-II	Buprenorphine
601622	VERDICT-II 7 Panel	AMP/BAR/BZO/COC/OPI3/PCP/THC

MEDTOX controls must be used with MEDTOX devices.