

VERDICT[®]-II DRUGS OF ABUSE TESTING

QUICK, ACCURATE RESULTS

Rely on VERDICT-II cassette devices for reliable drugs of abuse testing. Offered in 1- to 7-panel configurations, these devices deliver accurate screening results within minutes.

VERDICT-II Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485:2003
- Accurate screening results in minutes
- On-board adulteration strip available on some products
- Tests for up to 7 drugs of abuse in a single device
- Technical support services available for all MEDTOX Diagnostics, Inc. products

VERDICT-II Screening Thresholds

Amphetamine (AMP)	-	1000 ng/ml
Barbiturates (BAR)	-	200 ng/ml
Benzodiazepine (BZO)	-	300 ng/ml
Cocaine (COC)	-	300 ng/ml
Methamphetamine (MAMP)	-	1000 ng/ml
3,4-Methylenedioxy-N-		
Methylamphetamine (MDMA)	-	1500 ng/ml
Methadone (MTD)	-	300 ng/ml
Opiates (OPI3)	-	300 ng/ml
Oxycodone (OXY)	-	100 ng/ml
Phencyclidine (PCP)	-	25 ng/ml
Tricyclic Antidepressant (TCA)	-	300 ng/ml
Cannabinoids (THC)	-	50 ng/ml

VERDICT-II — Traditional Drug Cut-Off Level Cassette with Pipette

Part No.	Product Description	Panel Configurations
601584	VERDICT-II	THC
601585	VERDICT-II	Cocaine
601586	VERDICT-II	Opiates
601624	VERDICT-II	Benzodiazepines
601629	VERDICT-II	TCA
601598	VERDICT-II	Methamphetamines (MDMA)
601630	VERDICT-II	Methadone
601675	VERDICT-II	Oxycodone
601587	VERDICT-II 2 Panel	COC/THC
601593	VERDICT-II 4 Panel	AMP/COC/OPI3/THC
601645	VERDICT-II 4 Panel	COC/MAMP (MDMA)/OPI3/THC
601643	VERDICT-II 4 Panel	AMP/COC/OPI3/THC + LFAS
601594	VERDICT-II 5 Panel	AMP/COC/OPI3/PCP/THC
601638	VERDICT-II 5 Panel	AMP/BZO/COC/OPI3/THC
601640	VERDICT-II 5 Panel	AMP/COC/OPI3/MAMP (MDMA)/THC
601651	VERDICT-II 5 Panel	BZO/COC/MAMP (MDMA)/OPI3/THC
601639	VERDICT-II 5 Panel LFAS	AMP/COC/OPI3/MAMP (MDMA)/THC + LFAS
601622	VERDICT-II 7 Panel	AMP/BAR/BZO/COC/OPI3/PCP/THC

MEDTOX controls must be used with MEDTOX devices.