



EZ-SCREEN[®] (HI) & PROFILE[®]-IIA

DRUGS OF ABUSE TESTING

FLEXIBLE AND QUICK RESULTS.

With MEDTOX[®] drug testing devices, you have a choice of a cup or cassette. Both the EZ-SCREEN (HI) cup and the PROFILE-IIA cassette provide fast and accurate screening results in minutes.

EZ-SCREEN or PROFILE-IIA Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485
- Accurate screening results in minutes
- On-board adulteration strip (AS/LFAS) available on all EZ-SCREEN (HI) and PROFILE-IIA products
- Panel options include up to 11 tests
- Technical support services available for all MEDTOX Diagnostics, Inc. products

EZ-SCREEN and PROFILE-IIA Screening Thresholds

Amphetamine (AM/AMP)	- 1000 ng/ml
Barbiturates (BA/BAR)	- 200 ng/ml
Benzodiazepine (BZ/BZO)	- 300 ng/ml
Cocaine (CO/COC)	- 300 ng/ml
Methamphetamine (MA/MAMP)	- 1000 ng/ml
Methadone (MT/MTD)	- 300 ng/ml
Opiates (OP/OPI2)	- 2000 ng/ml
Oxycodone (OX/OXY)	- 100 ng/ml
Phencyclidine (PC/PCP)	- 25 ng/ml
Propoxyphene (PP/PPX)	- 300 ng/ml
Cannabinoids (TH/THC)	- 50 ng/ml

PROFILE-IIA Test Panel Options

Part No.	Product Description	Panel Configurations
601599	PROFILE-IIA 5 Panel	AMP/COC/OPI2/PCP/THC + LFAS
601607	PROFILE-II 7 Panel (BAR/BZO)	AMP/BAR/BZO/COC/OPI2/PCP/THC
601660	PROFILE-IIA 7 Panel	AMP/BAR/BZO/COC/OPI2/PCP/THC + LFAS
601673	PROFILE-IIA 7 Panel (MAMP/OXY)	AMP/COC/MAMP/OPI2/OXY/PCP/THC + LFAS
601677	PROFILE-IIA 10 Panel	AMP/BAR/BZO/COC/MAMP/OPI2/OXY/PCP/PPX/THC + LFAS

EZ-SCREEN Cup Test Panel Options

Part No.	Product Description	Panel Configurations
605211	EZ-SCREEN (HI) 4 Panel	AM/CO/OP/PC + AS
605200	EZ-SCREEN (HI) 5 Panel	AM/CO/OP/PC/TH + AS
605201	EZ-SCREEN (HI) 6 Panel	AM/CO/MA/OP/PC/TH + AS
605202	EZ-SCREEN (HI) 7 Panel	AM/CO/MA/OP/OX/PC/TH + AS
605203	EZ-SCREEN (HI) 8 Panel	AM/BA/BZ/CO/MA/OP/PC/TH + AS
605212	EZ-SCREEN (HI) 9 Panel	AM/BA/BZ/CO/MA/MT/OP/OX/PC + AS
605204	EZ-SCREEN (HI) 10 Panel	AM/BA/BZ/CO/MA/MT/OP/OX/PC/TH + AS
605205	EZ-SCREEN (HI) 11 Panel	AM/BA/BZ/CO/MA/MT/OP/OX/PC/PP/TH + AS

MEDTOX controls must be used with MEDTOX devices.