

PROFILE[®]-V for MEDTOXScan[®] DRUGS OF ABUSE TESTING

PERFORMANCE. EFFICIENCY. VALUE.

MEDTOX[®] PROFILE-V on-site devices bring greater clarity to testing for drug abuse. A reliable, fast and convenient alternative to in-lab testing, our MEDTOX PROFILE-V devices and MEDTOXScan reader provide accurate screening results.

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485
- Easy installation and operation
- One-step testing process
- Accurate screening results in ten minutes
- Lower clinical cutoff levels for AMP, BAR, BZO, COC, MAMP, MTD
- Eliminates human interpretation
- Load device into reader and walk away
- Technical support service available for all MEDTOX Diagnostics, Inc. products

MEDTOXScan Reader Features

- Stores up to 1,000 screening results with collector information
- Barcode scanner for accurate entry of patient identification, specimen ID and lot number
- Provides a digital display of screening results
- Printer provides paper copy of screening results
- Provides operator identifications
- Automated device detection
- Ability to interface with hospital laboratory information systems

Ordering Information for PROFILE-V Devices

Part No.	Product Description	Panel Configurations
604028	PROFILE-V 10 Panel	AMP, BAR, BZO, COC, MAMP, MTD, OPI2, OXY, PCP, THC

MEDTOX controls must be used with MEDTOX devices.



PROFILE-V Screening Thresholds

Amphetamine (AMP)	-	500 ng/ml
Barbiturates (BAR)	-	200 ng/ml
Benzodiazepine (BZO)	-	150 ng/ml
Cocaine (COC)	-	150 ng/ml
Methamphetamine (MAMP)	-	500 ng/ml
Methadone (MTD)	-	200 ng/ml
Opiates (OPI2)	-	2000 ng/ml
Oxycodone (OXY)	-	100 ng/ml
Phencyclidine (PCP)	-	25 ng/ml
Cannabinoids (THC)	-	50 ng/ml