



PROFILE-II ER[®]

DRUGS OF ABUSE TESTING

PERFORMANCE. EFFICIENCY. VALUE.

MEDTOX[®] PROFILE-II ER drug testing devices are accurate, easy to use and cost effective. Simply add the urine specimen to one of our cassette devices and get accurate screening results in minutes.

PROFILE-II ER Screening Thresholds

Amphetamine (AMP)	-	1000 ng/ml
Barbiturates (BAR)	-	200 ng/ml
Benzodiazepine (BZO)	-	300 ng/ml
Cocaine (COC)	-	300 ng/ml
Methamphetamine (MAMP)	-	1000 ng/ml
3,4-Methylenedioxy- Methamphetamine (MDMA)	-	1500 ng/ml
Methadone (MTD)	-	300 ng/ml
Opiates (OPI3)	-	300 ng/ml
Oxycodone (OXY)	-	100 ng/ml
Phencyclidine (PCP)	-	25 ng/ml
Propoxyphene (PPX)	-	300 ng/ml
Tricyclic Antidepressant (TCA)	-	300 ng/ml
Cannabinoids (THC)	-	50 ng/ml

PROFILE-II ER Features and Benefits

- 510(k)-cleared
- Made in the USA
- Quality management system conforms with ISO 13485
- Accurate screening results in minutes
- Panel options include up to 12 tests
- Technical support services available for all MEDTOX Diagnostics, Inc. products

Ordering Information for PROFILE-II ER Devices

Part No.	Product Description	Panel Configurations
601591	PROFILE-II ER 9 Panel	AMP, BAR, BZO, COC, MTD, OPI3, PCP, TCA, THC
601636	PROFILE-II ER 10 Panel	AMP, BAR, BZO, COC, MAMP-MDMA, MTD, OPI3, PPX, TCA, THC
601672	PROFILE-II ER 11 Panel	AMP, BAR, BZO, COC, MAMP-MDMA, MTD, OPI3, PCP, PPX, TCA, THC
601674	PROFILE-II ER 12 Panel	AMP, BAR, BZO, COC, MAMP-MDMA, MTD, OPI3, OXY, PCP, PPX, TCA, THC

MEDTOX controls must be used with MEDTOX devices.